

# ANALYSIS OF PERFLUOROALKYL SUBSTANCES (PFAS) AND METABOLITES IN BIOFLUIDS

**Waters**  
THE SCIENCE OF WHAT'S POSSIBLE.™

Mary Lane, *Kari Organtini*, *Beatriz Renner*, and *Ken Rosnack*  
Waters Corporation, Milford, MA USA

## INTRODUCTION

Perfluoroalkyl substances (PFAS) are common, man-made, persistent environmental contaminants that are used in the production of many consumer and industrial products. While highly advantageous for their chemical properties, these analytes don't metabolize or degrade easily, leading to increased persistence in the environment. Exposure and bioaccumulation of some PFAS residues has been linked to harmful health effects including: cancer, obesity, immune and reproductive issues.

With their high environmental prevalence, known bioaccumulation, and toxicity, these PFAS compounds have gained increased attention as global contaminants. Thus, environmental monitoring and biomonitoring of PFAS has steadily increased. With its many benefits (e.g., multiplexing, selectivity, dynamic range and fast method development), liquid chromatography-mass spectrometry (LC-MS) is the most common and widely accepted technique for their determination and highly sensitive and selective quantification.

However, developing robust methodologies for PFAS analysis can be quite challenging. The broad family of PFAS compounds have diverse chemical properties. Some can experience poor solubility (particularly of the longer chain PFAS) while others are strong anionic species with the propensity to adhere to various surfaces that they may come in contact with (e.g., sample containers, pipette tips, extraction plates/cartridges, collection vials/plates, and LC-components). Confounding these challenges are the low expected concentration levels in both serum and environmental matrices. Finally, known contamination from LC instrument components, solvents and common laboratory consumables makes LC and sample preparation method development difficult.

This work described herein provides a single, simple method for the simultaneous quantification of multiple PFAS substances from serum. This method uses analytical scale sub-2 $\mu$ m LC and fast, selective sample prep in 96-well format to achieve an recovery and low matrix effects, while achieving LODs as low as 0.5-1.0 ng/mL.

## METHODS

### Sample Preparation

Concentrated PFAS standards (legacy and emerging) were prepared as a mix in methanol and then spiked into fetal bovine serum (100  $\mu$ L). The pretreated sera samples were prepared for LC-MS/MS analysis using a Waters Oasis™ WAX  $\mu$ Elution 96-well SPE extraction plate and the protocol described below.

### Sample Extraction with Oasis WAX

**Condition:** 200  $\mu$ L 5% NH<sub>4</sub>OH in methanol followed by 200  $\mu$ L methanol

**Equilibrate:** 200  $\mu$ L water

**Load sample:** Entire diluted acid pretreated sample (~200  $\mu$ L) was loaded onto the extraction plate

**Wash 1:** 200  $\mu$ L acid in water

**Wash 2:** 200  $\mu$ L 0-10% acetonitrile or methanol in water

**Elute:** 2 X 25  $\mu$ L 1-5% NH<sub>4</sub>OH in 50-100% methanol or acetonitrile

**Dilute:** 50  $\mu$ L 2 mM ammonium acetate in water


**LC-MS analysis:** 10  $\mu$ L of prepared sample was injected for analysis.

### LC-MS Conditions


LC-MS/MS quantification of PFAS analogues was performed using a Waters Xevo™ TQ-S triple quadrupole MS (ESI-). Chromatographic separation was achieved using an ACQUITY™ H-Class UPLC PLUS system fitted with the PFAS Analysis Kit using an ACQUITY UPLC™ BEH C<sub>18</sub>, 2.1 x 100 mm, 1.7  $\mu$ m column (35 °C) at flow rate of 0.3 mL/min. Mobile phases used for analysis consisted of Water: Methanol (95:5) containing 2 mM ammonium acetate (Mobile Phase A) and Methanol containing 2 mM ammonium acetate (Mobile Phase B). Total analysis time was 22 minutes.

A total of 37 PFAS analogues were monitored for quantification. MRM conditions for each PFAS compound was optimized using MassLynx™ Software and the QuanOptimize tool. Full details on the LC/MS method, can be found in Waters Application Note 720006471EN.


## RESULTS



**Figure 1. Representative PFAS sample preparation method development and optimization using mixed-mode WAX SPE in the 96-well format.** Best recovery with minimal matrix effects (ME) was seen when using 1:1 pretreatment with 50% formic acid (FA) or 4% H<sub>3</sub>PO<sub>4</sub> in water (A). Washing the sample with a strong acid (10% FA or 4% H<sub>3</sub>PO<sub>4</sub>) maintained recovery and minimized MEs (B). MEs could be further reduced using a Wash 2 solution containing 2-10% organic without reducing recovery (C). Methanol containing 1-5% NH<sub>4</sub>OH provided best overall recovery of the PFAS analytes, while providing high selectivity with MEs < 15% (D).



**Figure 2. Representative matrix matched standard curves (0.5-20 ng/mL) for 6-2FTS (A), FHpPA (B), and GenX (C) extracted from 100  $\mu$ L bovine serum using mixed-mode WAX SPE in the 96-well format.**



**Figure 3. Representative LC-MS chromatograms, acquired using the ACQUITY H-Class UPLC and Xevo TQ-S Mass Spectrometer, for analysis of PFOS (A), PFOA (B), PFBA (C), and PFNA (D) extracted from 100  $\mu$ L of spiked sera (1-20 ng/mL) using mixed-mode WAX SPE in the 96-well format.**

## CONCLUSIONS

- An optimized SPE extraction method was developed for simultaneous extraction of PFAS analogs using a mixed-mode weak anion exchange polymeric sorbent, which afforded recovery and low matrix effects (Figure 1) ensuring sensitivity and selectivity of the assay for the diversity of polar and hydrophobic PFAS.
- Use of the  $\mu$ Elution SPE plate format eliminated the need for evaporation, reducing losses due to volatility, adsorption and re-solubilization. Sample preparation with simple SPE was < 30 minutes.
- A selective UPLC method employing a sub-2 $\mu$ m polymeric reversed-phase column, provided resolution from closely related PFAS compounds and endogenous interferences, improving overall selectivity and sensitivity.
- Use of a PFAS LC analysis kit (e.g., PEEK mobile phase tubing and PFAS isolator column) successfully delayed any potential PFAS contribution from the LC system.
- Detection or quantification limits of 0.5-1.0 ng/mL were achieved for various PFAS compounds extracted from 100  $\mu$ L of serum using a highly sensitive tandem quadrupole MS system.
- The method demonstrates it is fit-for-purpose for accurate measurement of PFAS compounds with high recovery and selectivity demonstrating its utility for monitoring exposure and potential health outcomes.